

To  
ACE Handels- und Entwicklungs GmbH  
Staufenstraße 1  
83395 Freilassing  
Deutschland / Germany

Herewith I/we give the order for the calibration of the following devices:

Model name	Amount	Calibration voucher included	
		yes	no
		yes	no
		yes	no
		yes	no
		yes	no

Order of accessories (in addition to the calibration order)		
Mouthpieces*	for breathalyzer model	Quantity
Replacement sensor	for breathalyzer model	Quantity

\*Standard PU = 25 pcs. / Dräger 3820: PU = 5 pcs.

**Address** ACE Handels- und Entwicklungs GmbH  
(ACE Instruments), Staufenstr. 1,  
D-83395 Freilassing / Germany  
**T** +49(0)8654 / 779 53 31 · **F** +49(0)8654 / 779 96 94  
**E-Mail** info@alkomat.net **Web** www.alkomat.net ·  
**CEO** Roland Heinemann · Klaus Forsthofer  
**Bank Account Germany**  
Postbank Munich · **Bank Code No.** 7001 0080 ·  
**IBAN** DE26 7001 0080 0804 3568 01 · **BIC** PBNKDEFF  
**Bank Account Austria**  
Salzburger Sparkasse · **Bank Code No.** 20404 ·  
**IBAN** AT982040406605002415 · **BIC** SBGSAT2S  
**VAT No. Germany** DE233838866 ·  
**Tax number Germany** 163/121/20473 ·  
**Commercial Register** HRB 18204  
(DISTRICT COURT Traunstein) ·  
**VAT No. Austria** ATU72268789 · **Tax number Austria** 430/9693  
**VAT No. Italy** IT00179419999 · **Spain** ESN2760709B ·  
**Great Britain** GB243240247 · **France** FR41820828010 ·  
**Poland** PL5263167208 · **Czech Republic** CZ683922218  
**www.alkomat.net** · **www.ace-technik.com**  
**www.gasmesstechnik.de** · **www.ace-brandschutz.com**  
**www.ace-markenshop.com** · **www.ace-arbeitsschutz.com**  
**www.msa-arbeitsschutz.com** · **www.ace-military.com**

## Calibration return receipt

<b>Invoice address</b>		ACE customer number (if available)
Company name		
First name	Last name	
Street / No.		
Postal code	City	
Country		
VAT number		
<b>Return address</b>		
Company name		
First name	Last name	
Street / No.		
Postal code	City	
Country		

<b>Contact data for call backs</b>	
Tel	
E-mail	
<b>Choose form of payment</b>	
<input type="checkbox"/>	Private individuals & Companies / Payment via PayPal**
<input type="checkbox"/>	Private individuals & Companies / Amount enclosed in cash
<input type="checkbox"/>	Private individuals & Companies / Payment via postal cash order
<input type="checkbox"/>	Companies / Payment via invoice (subject to examination)
**Please do the PayPal payment asap and inform us about the payment afterwards (by e-mail to info@alkomat.net or by phone call), thank you!	
Date / Company stamp / Signature of the client	

## Information about the calibration service

Print the calibration return receipt, fill in all required information and send it together with the breathalyzer (in a franked parcel) to the address to be seen on the Calibration return receipt. Please note that with the Calibration return receipt you have no claim for a post-free return.

This service is carried out at a fixed gross flat rate incl VAT and incl. return shipping in an insured DHL parcel.

Due to higher shipping costs, other flat rates apply for customers who do not live in AT or DE.

The respective flat rate price includes the following services and products for all customers (per device sent in):

- Calibration service with **return within 5 business days**
- Sensor inspection using the latest measuring instruments
- Professional calibration in the ACE service laboratory by technically top trained employees
- Attaching a sticker as reminder for the next calibration
- **Test report (calibration certificate)**
- **1 set of new batteries** (quantity depends on model) **and 25 additional mouthpieces** (only for mobile breathalyzers) **(exception: 5 mouthpieces for Dräger Alcotest 3820)**

## Why should a breathalyzer be regularly calibrated?

- In order to preserve the accuracy of breathalyzers the integrated sensor should regularly be readjusted = calibrated. The ultimate factor is not - as often assumed wrongly - the frequency of the use, but the time:
  - > For often used devices, the zero point is displaced. This means that after a breath test too high or too low values are displayed.
  - > In case of little used devices or a long period of time without using the device this also influences the sensitivity of the integrated sensor. It as well leads to too high or too low values after a breath test.
  - > So regardless of the use, a breathalyzer should be calibrated in certain time intervals (see details in the individual manual).

- **Guarantee expires if you don't pay attention to the calibration ranges.**

Manufacturer / Brand	Gross price for the calibration service	
	DE / AT	Other countries
Dräger Alcotest	40,- €	50,- €
EnviteC AlcoQuant	40,- €	50,- €
LION	40,- €	50,- €
ACE series (devices with fuel-cell sensor)	40,- €	50,- €
ACE series (devices with semiconductor sensor)	30,- €	40,- €
ACE series (stationary devices)***	80,- €	90,- €

\*\*\*For stationary devices (AL4000, ACE Public (V), ACE Vision, ACE Agent) the calibration costs are higher due to a more complex service.

For breathalyzers of the AL6000 and AL7000 series, for the devices ACE ALP-1 med and ACE III Premium as well as for the stationary breathalyzers you can purchase replacement sensors. ACE can also calibrate these replacement sensors. Should a repair or sensor renewal be necessary, we will contact you in advance. If you have any questions you may also contact our service staff by e-mail to info@alkomat.net or by telephone under +49 (0) 8654 77953-31.

Information on Revocation Instruction, Disclaimer, Data Protection Declaration, General Terms and Conditions see www.alkomat.net.

## Calibration at ACE Instruments

■ Each breathalyzer which you buy from ACE Instruments, can also be calibrated by ACE Instruments. Of course we also calibrate devices which were bought from other contractors, as long as these devices are identical in construction to the devices sold by ACE Instruments.

■ Most instruments give an indication of the recommended time for calibration. However, it can also be visible on a calibration sticker. If there is no calibration reminder and no calibration sticker our recommendation is that the calibration is done at latest 12 months after the purchase date and afterwards always 1x per year.

■ **Alcohol testers contain sensitive electronics: For this reason, please pack the device in such a way that no transport damage is possible. We recommend sending the device as a parcel so that it can be tracked if it cannot be found or that a proof of delivery can be provided.**